



SKU	Designation	French Law	Caliber	Weight (g)	Colors	MSRP
BB2027	Bag 1KG tracer green 0.28g ORGANIC BLS	Vente libre	6 mm airsoft	0.28	Green	35.30 € incl. tax
PCKBB2027	Box of 20 Bags 1KG trace green 0.28g ORGANIC BLS	Vente libre	6 mm airsoft	0.28	Green	599.00 € incl. tax

Enhance your nighttime games with high-precision, eco-friendly tracer balls.

BLS biodegradable tracer pellets are designed for demanding airsoft players looking to combine performance with environmental friendliness. Made from biodegradable PLA, they decompose naturally, reducing their environmental impact. Their 5.95mm diameter with a tolerance of ± 0.01 mm ensures smooth feeding and a stable trajectory, ideal for high-rate replicas. The green color of the pellets allows for optimal visibility during shots, making it easier to track the trajectory.

Technical characteristics:

- **Weight**
- **Caliber:**
- **Material**
- **Quantity**
- **Packaging**

Tracer Beads for an Unparalleled Nighttime Experience

BLS BIO Green Tracer BBs are specially designed for airsoft games in low light or night conditions. Thanks to their phosphorescent coating, they light up when fired, providing optimal visibility of the trajectory. This feature not only improves accuracy, but also adds a spectacular visual dimension to your games.

Eco-responsible Manufacturing

Made from **PLA (polylactic acid)** , a biodegradable polymer derived from natural resources such as corn starch, these beads decompose naturally without harming the environment. They are an eco-friendly alternative to traditional plastic beads, addressing growing concerns about sustainability and environmental impact.

Compatibility and Performance

With a precise diameter of **5.95 mm ± 0.01 mm** , these BBs ensure perfect compatibility with replicas equipped with tracer devices, such as tracer silencers or special hop-up chambers. Their weight of **0.20g to 0.28g** is ideal for replicas displaying an energy between 1.5 and 1.9 joules, offering a perfect balance between range and precision.

Directions for Use

To take full advantage of the tracer effect, it is recommended to use these balls with a tracer device that will activate their green luminescence when fired. If you do not have such a device, you can expose the balls to a light source before use to temporarily activate their phosphorescence. Store the balls in a dry place away from light to preserve their properties.

Why Choose BLS BIO Green Tracer Beads?

- **Superior Quality** : Impeccable finish and strict dimensional tolerance for optimum performance.
- **Respect for the Environment** : Made from biodegradable materials reducing the ecological impact.
- **Enhanced Gaming Experience** : Visible light path for immersive nighttime gameplay.
- **Wide Compatibility** : Suitable for a wide range of replicas equipped with tracer devices.

Les prix de vente conseillés sont mentionnés à titre indicatif. Les armuriers sont libres de vendre au prix qu'ils souhaitent. Textes et photos non contractuels, sujet à modification.